

vetOvation

## Skin Stapler



### Features and Benefits

The vetOvation Skin Stapler is a top quality, sterile, single use, skin stapler.

- **55 wide medical grade stainless steel staples** - most other brands have only 35 staples
- **Customer price \$12.00 each** - one stapler
- **Customer price \$120.00** - box of 10 staplers
- **Staple size: L 7.00mm, W 7.00mm, Dia: 060mm**
- **Unique ratchet system** - prevents jamming
- **Ergonomic, user friendly design**
- **Easy to check remaining staples in cartridge** - visible counter on side of instrument
- **Clear view on operating site**
- **Short trigger stroke** - for ease of use and nonstop feed.
- **Angled head allows visibility of incision** - ensures precise staple placement
- **Regular or wide staplers available**