


## Laparoscopic Instruments

Designed specifically for equine laparoscopy.

- 10mm shaft x 47cm long - with locking handle
- Babcock (non traumatic)
- Bowel grasper (non traumatic)
- Metzenbaum scissor
- Traumatic Ovary grasper
- 5mm shaft x 46cm long
- Similar instrument tips also available in a 5mm x 46cm version



## Vessel Sealer + Generator

The Aesculap Caiman is an easy-to-use vessel sealer that provides secure hemostasis for vessels up to 7mm in diameter during open and laparoscopic procedures.

- Distal tip instrument closure enables larger amounts of tissue captured inside the non-stick jaws
- Instruments average 50-75 procedures with EtO gas sterilization
- Shaft length 24, 36, 44cm
- Shaft diameter 5 & 12mm



## Trocars

Designed specifically for equine laparoscopy. Trocar cannulas come with magnetic valve and option of blunt or safety shield obturator.

- **Diameters: 5mm, 10mm, 12mm**
- **Shaft length: 21 cm**
- **Safety shield tip**
- **Sharp tip**
- **Blunt tip**
- **Visually guided**
- **Reducer cap available for 10mm and 12mm trocars**

Trocar with blunt obturator



Trocar with sharp obturator



Safety shield tip



Blunt tip



Reducer cap

