



Exclusive

PATTERSONVET.COM | 800.279.6452

**INTRODUCING
THE VETO VATION
VIDEO SCOPE.**
It works with your iPhone.





PATIENT CARE IN THE PALM OF YOUR HAND.

Patterson's exclusive vetOvation video scope is an affordable option for otology and rhinoscopy. It connects to your iPhone and uses its camera to record videos and still images to pinpoint problems inside the ear and nose. It's easy to share these videos and images with clients, too, so you can educate them by actually *showing* them the problem.



**ADAPTS TO YOUR
IPHONE OR IPOD TOUCH!**

iPhone® is a trademark of Apple Inc., registered in the U.S. and other countries.

BETTER PATIENT CARE THROUGH INNOVATION.

320 Lumen Light with
recharging station
rechargeable lithium battery

Adapter for iPhone 4s, 5/5s,
6/6+, 7 and iPod touch



Flexible serrated
grasper, 5 Fr x
260 mm

Scope 0 degrees
5.0 mm x 90 mm
Autoclavable,
dual channel
sheath accepts
serrated grasper
and continuous
lavage

THIS EASY-TO-USE DEVICE ENABLES:

- Video otoscopy for quick and easy inspection of the middle ear, removal of excessive hair and wax, along with examination of the tympanic membrane
- Rhinoscopy; the 5.0 mm 0-degree scope enables removal of foreign objects and biopsies
- Intubation of rabbits and other small animals
- Removal of foxtails

**Download the vetOvation
video scope app free at
the App Store.**

