

INTRODUCING THE VETOVATION SURGICAL SYSTEM.

*Minimally invasive
procedures are
now affordable.*



Exclusive

PATTERSONVET.COM | 800.279.6452



AN AFFORDABLE WAY TO EXPAND PATIENT CARE.

Patterson's exclusive vetOvation Surgical System is an affordable way for your practice to start offering minimally invasive procedures, including laparoscopic spay, prophylactic gastropexy, liver or kidney biopsy, cystoscopy, cryptorchid resection, rhinoscopy, cholecystectomy, otoscopy, soft-pallet resection and more.

Expand your practice with minimally invasive surgery options. The outcome? Shorter surgery and anesthesia times, healthier patients and happier clients.



OUR SURGICAL SYSTEM INCLUDES THE FOLLOWING:

- High-definition camera and monitor, CO2 insufflator and LED light, which are manufacturer refurbished from a human hospital. They come with a one-year warranty and are in excellent working condition.
- Tilt table attaches to any operating room table and can be rotated side to side.
- Vessel sealing technology for secure hemostasis for vessels up to 7 mm in diameter during open and laparoscopic procedures. Reduce the duration of anesthesia and decrease blood loss and tissue trauma.
- Laparoscopic instruments designed specifically for veterinary laparoscopy are available in two lengths and are autoclavable.
- Trocars are 5 mm in diameter. They're low profile with stability threads, making them perfect for small and large animals. Fully autoclavable.
- Two days of in-hospital training are also included with your purchase.



BETTER PATIENT CARE THROUGH EASY INNOVATION.