


Laparoscopic Instruments

Designed specifically for equine laparoscopy.

- 10mm shaft x 47cm long - with locking handle
- Babcock (non-traumatic)
- Bowel grasper (non-traumatic)
- Metzenbaum scissor
- Traumatic Ovary grasper
- Similar instrument tips also available in a 5mm x 46cm version
- 5mm x 35cm Equine Injection Needle



Vessel Sealer + Generator

The Voyant Intelligent Energy System by Applied Medical is an easy-to-use vessel sealer that provides secure hemostasis for vessels up to 7mm in diameter during open and laparoscopic procedures.

- Cuts and coagulates tissue during open and laparoscopic procedures using radiofrequency energy
- Instruments average 50-75 procedures with EtO gas sterilization
- Starter Kit includes generator and three Fusion instruments



Trocars

Designed specifically for equine laparoscopy. Trocar cannula come with a magnetic valve and the option of a safety shield, sharp or blunt obturator.

- **Diameters:** 5mm, 10mm, 12mm & 13mm
- **Shaft length:** 21 cm
- **Sharp tip**
- **Blunt tip**
- **Safety shield tip**
- Reducer cap available for 10mm, 12mm & 13mm trocars

Trocar with blunt obturator



Trocar with safety shield obturator



Safety shield tip



Blunt tip



Reducer cap

