




# Alexis<sup>®</sup>

## Wound Protector-Retractor


# PROTECT EVERY INCISION WITH ALEXIS WOUND PROTECTORS

## Alexis wound protector-retractors:

- Provide 360° of wound protection and atraumatic retraction.
- Maintain moisture at the incision site.
- Distribute force evenly, eliminating point trauma and associated pain.
- Free up valuable hands in the operating room.
- Accommodate incision sizes of 1-25cm.



Foreign object removal performed in a canine with Alexis wound protector-retractor.



Splenectomy procedure performed in a canine with Alexis wound protector-retractor.



Alexis O Wound Protector-Retractor

Alexis Wound Protector-Retractor





## Alexis O Wound Protector-Retractor

*Featuring a rigid retraction ring for maximum exposure*

Reorder No.	Size	Sheath Length	Incision Range
C8401*	Small	18cm	2.5-6cm
C8402	Medium	18cm	5-9cm
C8403	Large	25cm	9-14cm
C8404	X-Large	34cm	11-17cm
C8405	XX-Large	36cm	17-25cm



## Alexis Wound Protector-Retractor

*Featuring a flexible retraction ring for anatomical conformity*

Reorder No.	Size	Sheath Length	Incision Range
C8313*	XX-Small	20cm	1-3cm
C8323*	XX-Small, Short	11cm	1-3cm
C8312*	X-Small	19cm	2-4cm
C8322*	X-Small, Short	13cm	2-4cm
C8301*	Small	18cm	2.5-6cm
C8302	Medium	18cm	5-9cm
C8303	Large	25cm	9-14cm
C8304	X-Large	34cm	11-17cm

\*Model includes a tether to facilitate device removal.

Visit [www.vetovation.com](http://www.vetovation.com) for more information.

© 2019 Applied Medical Resources Corporation. All rights reserved.  
 Applied Medical, the Applied Medical logo design and marks designated with a ® are trademarks of Applied Medical Resources Corporation, registered in one or more of the following countries: Australia, Canada, Japan, South Korea, the United States, and/or the European Union.  
 VetOvation Affordable Minimally Invasive Solutions is a trademark of VetOvation, Inc.  
 251365-EN-USA-A

**vetovation**  
 Affordable Minimally Invasive Solutions



**ANIMAL HEALTH**