



Innovative
Products Designed
To Deliver
**Exceptional
Outcomes**

Unleash Your Precision

INTRODUCING


Hall® MicroFree®

Cordless Small Bone Power System

Designed to deliver comfort and simplicity with the performance and reliability you expect from Hall®.

- Cordless handpiece is designed to improve hand positioning and maneuverability.
- Well-balanced, ergonomic handpieces are designed to perform comfortably in surgeons' hands.²
- Adjustable pencil grip designed for precise cutting.

CONMED.com/MicroFree





HALL[®]
POWERED INSTRUMENTS

Control and Precision in the Palm of Your Hand



Hall® MicroPower+ Small Bone Handpiece System

Designed to deliver comfort and simplicity with the performance and reliability you expect from Hall®.



- Encapsulation and potting protect against moisture.
- Rotating sagittal saw head for improved cutting angles.
- New extendable levers with padded tip for in-hand comfort.
- Flared bur guards with silicon padding for handling and stabilization.
- Indexing oscillating saw collet design for cutting versatility.



Drill

Oscillating
Saw

Reciprocating
Saw



HALL[®]
POWERED INSTRUMENTS

Dependable Power that Works the Way You Do

Hall® 50™

Powered Instrument System

Designed to deliver comfort and simplicity with the performance and reliability you expect from Hall®.

- Hall's UL-approved autoclavable small and large lithium batteries deliver greater torque and longer performance than similarly sized batteries.*
- Handpieces are manufactured to provide a seamless and controlled clinical experience. They have also been engineered to produce less vibration, giving you better command of the instrument.
- A Hall 50 handpieces have been ergonomically designed to maximize surgeon comfort. They fit comfortably in your hands and are also light-weight which helps minimize fatigue during lengthy usage.
- Made in the USA and built by hand, every single handpiece goes through an exhaustive and comprehensive inspection by eight separate quality control experts before it ever leaves our factory.



HALL®
POWERED INSTRUMENTS

* Data on File.

▶ The Simplest High-Performance Pump in Arthroscopy

10k[®] Pump

High-Performance Inflow Only Pump

- Ease of use and superior performance are the key characteristics that make the 10k Fluid System unique.
- Quick setup and priming with intuitive one-step/one-hand cassette loading.
- Integrated pole mount works with any standard I.V. pole and compatible with video cart shelf or any flat surface.



Deliver new levels of simplicity, versatility and performance

D4000[®] Arthroscopic Resection System

- A versatile console that supports 31 handpieces and 5 languages. Features built-in wireless connectivity to support multiple wireless foot controllers and allows for customizable user-profiles that let users save their preferred settings.
- Designed to provide intuitive functionality with dual touchscreen interface so users can quickly and simply adjust a wide range of settings.
- Directly integrates with ConMed's 24k[®] Pump. These two units will effortlessly work together to automatically detect and adjust fluid levels.





Ergo™ Shaver Handpiece

Ergonomically engineered handle which provides optimal balance and control when held traditionally (knee), dagger style (shoulder and hip) or pencil grip. This eliminates the fatigue associated with traditional barrel shaped styles.

- Features: Bi-Directional Blade Insertion
- Weight: 19.3 oz. evenly distributed
- Max Speed: 12,000 RPM
- Torque: 35 in-ounces



Ultra Series™ Dual-Purpose Blades

- Can resect tissue and debride bone
- Feature a stronger, more robust inner cutting tip material and patented XtraSharp™ technology
- Superior cutting performance, longer region of lubrication and optimized dimensions



Better Outcomes Start With Better Visualization

INTRODUCING

4K Platform Solutions

- Sony 4K Monitor
- Ultra High-Definition Scopes
- CMOS Camera Heads, for greater visibility
- CMS9000™ Content Management System



CMS9000™ Content Management System

- Customizable surgeon settings
- Live streaming video
- Wireless connectivity and integration to hospital's EHR

Wireless tablet sold separately.



VP8500 Integrated Video Cart

- Compact and mobile, the VP8500 helps you use valuable space more efficiently.
- Five adjustable 14" shelves
- Fold-up keyboard tray
- Adjustable height 10k pump side mount
- Footswitch bin for extra storage



HALL[®]
POWERED INSTRUMENTS

mtfbiologics

 **Linvatec**

 **CONCEPT**[®]

Shutt[®]

CONMED Corporation
525 French Road, Utica, NY 13502

Local 727-392-6464
Toll Free 1-800-237-0169

CONMED.com
customer_service@CONMED.com