





# **Voyant**Intelligent Energy System

The Voyant advanced bipolar system collects information about the nature of the tissue within its jaws, rapidly and constantly measures tissue as the energy is applied, and adjusts to provide the optimal amount of energy throughout the seal cycle to create a permanent, fused seal.



Leg amputation performed in a canine using the Voyant Fine Fusion device.



Splenectomy procedure performed in a canine with the Voyant Open Fusion device.



# Device Key with Embedded Intelligence

Stores activation data to learn from live tissue.

Optimizes energy delivery for procedural and patient needs.

Enables efficient implementation of algorithm updates by delivering the most advanced technology with new handpieces.

# **Voyant Maryland Fusion**



#### **Jaw Features**

20mm seal length 18mm cut length

Single-action, curved and tapered jaw with dissecting tip

Single-handed 360° shaft rotation



# **Voyant 5mm Fusion**



#### **Jaw Features**

20mm seal length 18mm cut length Single-action, straight jaw with blunt tip Single-handed 360° shaft rotation



# **Voyant Open Fusion**



#### **Jaw Features**

40mm seal length
38mm cut length
Dual-action, curved jaw with blunt tip
180° shaft rotation



# **Voyant Fine Fusion**



#### **Jaw Features**

17mm seal length
15mm cut length
Dual-action, curved jaw with dissecting tip





## **Laparoscopic Devices**

Reorder No.	Description	Shaft Dimensions
EB210	Voyant 5mm Fusion	5mm x 37cm
EB211	Voyant 5mm Fusion	5mm x 44cm
EB215	Voyant Maryland Fusion	5mm x 37cm
EB216	Voyant Maryland Fusion	5mm x 44cm
EB217	Voyant Maryland Fusion	5mm x 23cm



## **Open Devices**

Reorder No.	Description	Shaft/Device Length
EB230	Voyant Fine Fusion	19.3cm
EB240	Voyant Open Fusion	20cm



#### Generator

Reorder No.	Description
EA020	Voyant electrosurgical generator

Visit www.vetovation.com for more information.



vet@vation