


# Vessel Sealing Technology

Reduce Tissue Trauma, Blood Loss and Surgery Time



OPEN FUSION

FINE FUSION

5MM FUSION

**PROCEDURES:** SPAY | NEUTER | PYOMETRA RESECTION | SPLENECTOMY | LEG AMPUTATION  
MAMMARY CHAIN MASS REMOVAL | SOFT PALATE RESECTION

## Vessel Sealer

The Voyant vessel sealer is designed for a variety of open and laparoscopic procedures, providing exceptional coagulation and transection of vessels up to 7mm in diameter.

- Voyant Electrosurgical Generator
- Multiple Voyant Vessel Sealing Devices
  - Voyant Fine Fusion
  - Voyant Open Fusion
  - Voyant 5mm Fusion
  - Voyant Maryland Fusion

*Low-temperature sterilization required.*

ITEM	VETOVAATION CODE	PATTERSON CODE	PATTERSON PRICE
Voyant Generator	EA020	78937903	\$11,500.00
Open Fusion Instrument	EB240	78937905	\$720.00
Fine Fusion Instrument	EB230	78937904	\$720.00
5mm Fusion Instrument	EB210	78937906	\$720.00
Maryland Fusion 23cm	EB217	78940880	\$720.00



# Sterlink Sterilizer

Low-temperature hydrogen peroxide sterilization allows veterinary practices the ability to sterilize many single patient use items along with batteries, drills, vessel sealing handles, scopes and camera heads. Sterilization time from 7–36 minutes.

## Low-Temperature Hydrogen Peroxide Sterilizer

- Safe low-temperature alternative to ethylene oxide
- Cycle time only 7, 14 or 36 minutes!
- Non-toxic and eco-friendly

### STERLINK Sterilizer + Supply Starter Kit



ITEM	VETOATION CODE	PATTERSON CODE	PATTERSON PRICE
STERLOAD® Full Chamber Cartridge (Box of 30)	STER-LOAD	78938639	\$240.00
STERPACK Plus® Medium Pouch (Box of 30)	STER-PLUS	78938645	\$210.00
STERPACK® Small Pouch (Box of 50)	STER-PACK	78938641	\$300.00
STERLINK® FPS + Starter Kit	STER-KIT	78938928	\$14,500.00



# Electrosurgery & Smoke Evacuation

Specialized veterinary radio frequency electrosurgical systems that communicate automatically with a smoke evacuation and filtration unit for safety and ease of use.

## Electrosurgery Kit includes:

- Reusable monopolar pencil
- Electrode tips
- Foot pedal
- Reusable grounding electrode/cable

ITEM	VETOATION CODE	PATTERSON CODE	PATTERSON PRICE
Surtron 200 Electrosurgery Generator	VES-ESU-00435	78947675	\$5,000.00
Surtron 160 Electrosurgery Generator	VES-ESU-00436	78947676	\$1,800.00
Surtron Smoke Evacuator	VES-ESU-00442	78947682	\$1,900.00



**200 Watt**— Digital high frequency electrosurgical generator that performs precision monopolar and bipolar surgery in addition to microsurgery. It is designed for ease of use utilizing a touch screen interface. This electrosurgery unit offers the ability for hand or foot control.



**160 Watt**— Monopolar radio frequency electrosurgical generator that performs CUT, Forced or Soft COAG for superficial or deep penetration and BLEND for a combination of cut and coagulation. There is also an option for BIPOLAR coagulation with a special adapter. This electrosurgery unit offers the ability for hand or foot control.



**Smoke Evacuation**— An aspiration and filtration system for the smoke produced during a surgical procedure. The system communicates seamlessly with the 160 & 200 watt electrosurgery units for safety and ease of use.