

vetOvation Video Otoscope



An affordable scope for video otology and rhinoscopy. VetOvation's video otoscope connects to an iPhone™, using its camera technology to process videos and still images to pinpoint recurrent or underlying problems inside the ear and nose.

This easy-to-use device enables:

- **Video Otoscopy** for animals with recurrent ear infections. The system allows for quick and easy inspection of the middle ear, removal of excessive hair and wax, and examination of the tympanic membrane.
- **Rhinoscopy** - the 6.0mm 0o scope enables removal of foreign objects and biopsies.
- The vetOvation video otoscope also allows sharing of videos and images with clients through email.



Components

- **Scope** - 0° 6.0mm x 90mm autoclavable; dual channel sheath accepts serrated grasper and continuous lavage
- **Flexible serrated grasper** - 5Fr x 260mm
- **320 Lumen light** with rechargeable lithium battery
- Adapters available for iPhones™ : 10, 10+, 7/7s & 8/8s, 7+, 6/6s, 6+, 5/5s, 4s, iPod Touch
- Free iPhone™ app



"We all have smart phones. Now with the VetOvation video otoscope we are able to view anatomy inside areas our eyes just can't see. Having the ability to capture real time pictures and video enables better treatment for the patient and increased compliance by the pet owner."

*-Dr. Chuck Miller Practice Owner
Triangle Veterinary Hospital, Durham, NC*