

IntelliScope WiFi HD Camera System




HD WiFi wireless camera system attaches to any standard scope. Share the image to your smart phone, tablet or PC.

- Video Otoscopy
- Rhinoscopy
- Cystoscopy
- Arthroscopy
- Intubation





WiFi Camera Handle
\$1,500
Patterson Product Code: 78947659
Vetovation Product Code: VID-WIF-00429



Video Otoscopy System

- **Scope** – 0 degree 6.0mm x 90mm, dual channel sheath accepts grasper + lavage
- **320 Lumen LED light** with charging base
- **Flexible serrated grasper** – 5Fr x 250mm
- **HD WiFi camera head** with coupler

\$2,875

Patterson Product Code: 78947670
Vetovation Product Code: VID-WIF-00430



Interchangeable Otoscope/ Rhinoscope System

- **Scope** – 0 degree 2.7mm x 110mm
- **Otoscope cone** with dual channel adaptor
- **Rhinoscope sheath**
- **320 Lumen LED light** with charging base
- **Flexible serrated grasper** – 5Fr x 250mm
- **HD WiFi camera head** with coupler

\$3,375

Patterson Product Code: 78947671
Vetovation Product Code: VID-WIF-00431



Multipurpose Scope System

- **Scope** – 30 degree 2.7mm x 188mm
- **Dual channel sheath** accepts grasper + lavage
- **320 Lumen LED light** with charging base
- **Flexible serrated grasper** – 5Fr x 250mm
- **HD WiFi camera head** with coupler

\$3,875

Patterson Product Code: 78947672
Vetovation Product Code: LAP-WIF-00432

TCI Scope System

\$7,250

Vetovation Product Code: VID-WIF-00988

Arthroscopy Kit-1.9mm Scope

\$4,275

Patterson Product Code: 78947673
Vetovation Product Code: ORT-WIF-00433

Arthroscopy Kit-2.5mm Scope

\$4,275

Patterson Product Code: 78947674
Vetovation Product Code: ORT-WIF-00434

vetOvation Video Otoscope



An affordable scope for video otology and rhinoscopy. VetOvation's video otoscope connects to an iPhone™, using its camera technology to process videos and still images to pinpoint recurrent or underlying problems inside the ear and nose.

This easy-to-use device enables:

- **Video Otoscopy** for animals with recurrent ear infections. The system allows for quick and easy inspection of the middle ear, removal of excessive hair and wax, and examination of the tympanic membrane.
- **Rhinoscopy** - the 6.0mm 0o scope enables removal of foreign objects and biopsies.
- The vetOvation video otoscope also allows sharing of videos and images with clients through email.



Components

- **Scope** - 0° 6.0mm x 90mm autoclavable; dual channel sheath accepts serrated grasper and continuous lavage
- **Flexible serrated grasper** - 5Fr x 260mm
- **320 Lumen light** with rechargeable lithium battery
- Adapters available for iPhones™ : 13 Pro Max, 13/13 Pro, 12 Pro Max, 12/12 Pro, 11 Pro Max, 11/11 Pro, XR, XS, 10, 10+, 7/7s & 8/8s, 7+, 6/6s, 6+
- Free iPhone™ app

Video Otoscope Kit

\$1,750



"We all have smart phones. Now with the VetOvation video otoscope we are able to view anatomy inside areas our eyes just can't see. Having the ability to capture real time pictures and video enables better treatment for the patient and increased compliance by the pet owner."

*-Dr. Chuck Miller Practice Owner
Triangle Veterinary Hospital, Durham, NC*

