

STERLINK®

STERLINK U-510



VETERINARY EXCELLENCE REDEFINED

Stay at the forefront of veterinary care with Plasmapp's cutting-edge tabletop plasma technology, providing a solution that aligns with the latest advancements in the field



STERLINK U-510

- ACHIEVE COMPLETE STERILIZATION CYCLES IN JUST 17 OR 24 MINUTES
- MAXIMIZE STERILIZATION CAPABILITIES WITH AN IMPRESSIVE 20 CYCLES
- MINIMIZE PRODUCT USAGE WHILE MAINTAINING THE HIGHEST STANDARDS OF STERILIZATION

STERLINK

STERLINK U-510

TECHNICAL DATA SHEET



SPECIFICATIONS

STERLINK U-510

Total size (W x D x H)	21.34 x 16.46 x 11.93 inches
Weight	66.14 pounds
Capacity	2.06 gallons
Chamber (W x D x H)	18.5 x 7.28 x 3.54 inches
Diagonal Length (Chamber)	20.1 inches
Vacuum Pump	Pump built-in type
Cycle Time	Standard Mode: 17 minutes Advanced Mode: 24 minutes

VAPORIZED STERILANT CIRCULATING SYSTEM

The star of the U-510 experience is the STERLOAD, a remarkable accessory that extends your sterilization capabilities. Its leak-proof double construction releases gases from the natural decomposition of the sterilant. With the STERLOAD, you can now achieve an incredible 20 cycles in just 17 or 24 minutes, maximizing efficiency and streamlining your workflow.

