


360° Wound Protection

Makes Wound Protection and Retraction Hands-Free



The Alexis Wound Protector-Retractor:

- 360° Wound Protection
- Reduces Surgical Site Infections (SSIs)
- Atraumatic, Hands-Free Retraction
- Better Visualization
- Minimizes Blood Loss



Alexis

- Only option for patients requiring the XX-Small or X-Small sizes
- Preferred when the flexible ring's ability to conform to anatomy is an advantage
- For irregular or curved anatomy or thin body walls
- Conformity and gentle retraction
- Used more for small animals

Alexis O Rigid

- Only option for patients requiring the XX-Large size
- Preferred when maximum rigid exposure of a larger surgical field is needed
- For thick body walls and deep cavities
- Maximum rigid opening for deep visualization
- Used more in large animals

The Alexis Wound Retractor



ORDER CODE	SIZE	PRICE	INCISION RANGE
C8313	XX-Small	\$60.00	1-3 cm
C8312	X-Small	\$60.00	2-4 cm
C8301	Small	\$60.00	2.5-6 cm
C8302	Medium	\$80.00	5-9 cm
C8303	Large	\$80.00	9-14 cm
C8304	X-Large	\$95.00	11-17 cm

Procedures: Exploratory Laparotomy | Foreign Body Removal | Cystotomy | Exotic Animal
Abdominal Surgery | Spay | Thyroidectomy | Pyometra Resection |
Splenectomy

The Alexis O Rigid Wound Retractor



ORDER CODE	SIZE	PRICE	INCISION RANGE
C8401	Small	\$60.00	2.5-6 cm
C8402	Medium	\$80.00	5-9 cm
C8403	Large	\$95.00	9-14 cm
C8404	X-Large	\$95.00	11-17 cm
C8405	XX-Large	\$95.00	17-25 cm

Procedures: Cesarean Sections | Equine Standing Flank Laparotomy | Large Animal
Exploratory Laparotomy | Large Animal Splenectomy | Canine Intestinal
Resection & Anastomosis | Canine Gastrotomy