

LeiLOX 1.5 / 2.0 Micro




1.5 Osteosynthesis Locking Set, Titanium
142-1500-00




2.0 Osteosynthesis Locking Set, Titanium
142-2000-00

1.5 / 2.0 COMBI SETS
Also available in our shop




1.5 / 2.0 Osteosynthesis Locking Set
142-1520-00

1.5 / 2.0 Challenger Locking Set
142-1520-50



TTA RAPID® Premium Set
132-6002-00



TTA RAPID® Starter Set II
132-6004-00


OTHER TTA
RAPID®SETS
Available in our shop




TTA RAPID® Starter Set III
132-6005-00

TTA RAPID® Starter Set IV
132-6006-00

LeiLOX 2.0 / 3.5



2.0 / 2.4 Straight Implant Set
142-2024-00



2.7 / 3.5 Reconstruction Implant Set
142-3100-00


TTA RAPID® Systems

Tiny TTA RAPID®



Tiny TTA RAPID® Set
132-6500-00

RAPID LUXATION



RAPID LUXATION Set
132-8000-01

High biocompatibility. Improved shape.

TPLO SWING IMPLANTS

A new product from RITA LEIBINGER & vetOvation



Now in
TITANIUM



Visit our online shop
Express shipping available

www.vetOvation.com


Or scan the QR-Code to
SHOP ONLINE




SHOP NOW

vetOvation

LeiLOX TPLO Swing



2.0 / 2.4 TPLO Swing Set
142-2820-24




2.7 / 3.5 TPLO Swing Set
142-2827-35

LeiLOX TPLO Locking



2.0 / 2.4 TPLO Set
142-0150-00



2.7 / 3.5 TPLO Set
142-0100-02

FEATURES & BENEFITS

Improved anatomic shape

Contoured to match the anatomic shape of the bone for better fit and positioning, and secure fixation

Very tight compression

Precisely designed double compression holes for a very tight compression of the osteotomy

Versatility through interchangeable screws

With identical screw heads, all Titanium 2.7/3.5 LeiLOX plates (TPLO and fracture systems) work with 2.7mm as well as 3.5mm screws in all of the plate sizes. Same applies for the 2.0/2.4 systems.



Titanium for best biocompatibility and durability

Made of medical grade titanium, which makes them very durable but lightweight, highly biocompatible, and low in temperature sensitivity

Multiaxial Locking for better flexibility

The screws can be locked in a 90° angle with a 12° deviation in any direction. This allows you to angle the screws away from vital structures.

LeiLOX CBLO Locking



2.0 / 2.4 CBLO Set
142-4200-00



2.7 / 3.5 CBLO Set
142-4000-10

Learn from the EXPERTS

Learn from the best and hone your skills in one of our Continuing Education (CE) Courses / Workshops, Online Courses, and Webinars.

Our **CE Courses and Workshops** are held by renowned experts and clinicians at locations all over the world, and cover everything from the fundamentals to advanced procedures.

Take a look at the range of educational courses we offer at the RITA LEIBINGER Academy website page.

Visit our Academy Page
leibinger.vet/academy

Contact us
contact@vetovation.com



ACADEMY

ALSO NEW

FROM RITA LEIBINGER MEDICAL

Suture Anchors

Made of Titanium and designed with self-tapping threaded tips for ease of use and a very secure and stable fixation to the bone.



LeiCOM

Our new range of Cannulated Compression Bone Screws in Titanium is available in Micro, Mini, Medi and Maxi sizes, and are also sold in sets.

vetovation

To contact a representative or receive additional information, please contact VetOvation.
919.247.0328 ■ contact@vetovation.com ■ vetovation.com